

REVLINE®

HOPPER THROTTLE™ DRY BEAN PLANTER BOX TREATMENT



NUTRIENT & MICROBIAL PLANTER BOX DELIVERY SYSTEM

Revline® Hopper Throttle® Dry Bean is an advanced planter box system that combines seed fluency, micronutrients and live biology in a convenient delivery system that can replace any seed fluency agent.

The patented Bio-Capsule™ Technology delivery system allows for the addition of multiple biological solutions safely packaged for convenient, live deployment at planting.

Treats the following: Black beans, Navy beans, Great Northern beans, Snap beans, Scarlet Runner beans, Wax beans, Pinto beans, Pink beans, Kidney beans, Garden or String beans, Field or Canning beans, Cranberry beans

PRODUCT BENEFITS

- **IonLock™ Zinc increases singulation & adheres microbes to seed coats while providing valuable early emergence nutrition**
- **Increases early-season iron, manganese and zinc uptake**
- **Drives early-season taproot nodulation and N-fixation**
- **Builds stronger, healthier plants fast & massive roots**
- **Puts dry beans in the position to manage season-long stressors and optimize genetic potential**


BioFertility & N-Fixing Microbes

Industry-Leading Rhizobium Leguminosarum Biovar Phaseoli Inoculant
 Accelerates early nodulation

Revline
 Proven, industry-leading biostimulant PPFM strains that generate massive root structures

Base Pail & Seed Fluency:

- 80/20 Talc/Graphite
- Active Carbon
- Mn & Fe
- 0.38 lbs IonLock Zinc equiv. to 1.35 qts 9% Zn EDTA



Actual product packaging may be different than this visual

GUARANTEED ANALYSIS

BASE INGREDIENTS

Iron (Fe).....	0.70%
Manganese (Mn)	0.90%
Zinc (Zn).....	1.58%
1.58% Water Soluble Zinc	
<i>Derived from Iron Oxide, Manganese Oxide and Zinc Sulfate.</i>	

BIO-CAPSULE INGREDIENTS:

ALSO CONTAINS BENEFICIAL SUBSTANCES

<i>Azospirillum brasilense</i>	5.0x10 ⁷ CFU/g
<i>Azospirillum lipoferum</i>	5.0x10 ⁷ CFU/g
<i>Azotobacter chroococcum</i>	5.0x10 ⁷ CFU/g
<i>Azotobacter vinelandii</i>	5.0x10 ⁷ CFU/g
<i>Bacillus amyloliquefaciens</i>	5.3x10 ⁸ CFU/g
<i>Bacillus licheniformis</i>	1.1x10 ⁹ CFU/g
<i>Bacillus megaterium</i>	5.3x10 ⁸ CFU/g
<i>Bacillus pumilus</i>	8.8x10 ⁸ CFU/g
<i>Bacillus subtilis</i>	5.3x10 ⁸ CFU/g
<i>Paenibacillus azotofixans</i>	5.0x10 ⁷ CFU/g
<i>Thiobacillus ferrooxidans</i>	5.0x10 ⁷ CFU/g
<i>Trichoderma harzianum</i>	3.5x10 ⁵ CFU/g
<i>Methylobacterium hispanicum</i>	1x10 ⁹ CFU/g
<i>Rhizobium Leguminosarum Biovar Phaseoli</i>	1.0x10 ⁵ CFU/g

DIRECTIONS FOR USE & RECOMMENDED USE RATE

DO NOT ACTIVATE BIO-CAPSULE UNTIL READY FOR USE. Read and follow all label directions for use on package. Once ready to plant, remove the safety clips from the Bio-Capsules. Push down on the buttons atop the pail to release the biologicals into the base. Seal the pail and shake aggressively to blend the contents. Before each application, aggressively shake contents to ensure the correct volume is applied.

The combined products (1 pail) treat 40 units of dry beans.

For single-row unit planters, use the enclosed scoop to measure out blended contents. One scoop treats 1 unit dry beans.

PRODUCT DETAILS

Sold in 2-pail cases to treat 80 units of dry beans. Each pail is capable of treating 40 units of dry beans.

©2025 Meristem Crop Performance® Group, Inc. All rights reserved. Bio-Capsule, Hopper Throttle, IonLock, and Revline are trademarks of Meristem Crop Performance Group, Inc. v20250929.

