

REVLINE®

HOPPER THROTTLE™ SOYBEAN PLANTER BOX TREATMENT



NUTRIENT & MICROBIAL PLANTER BOX DELIVERY SYSTEM

Revline® Hopper Throttle® Soybean Ether™ is an advanced planter box system that combines seed fluency, micronutrients and live biology in a convenient delivery system that can replace any seed fluency agent.

The patented Bio-Capsule™ Technology delivery system allows for the addition of multiple biological solutions safely packaged for convenient, live deployment at planting.

Ether Enzyme Technology kick-starts the uptake of available soil nutrients and water while maximizing microbial activity.

PRODUCT BENEFITS

- **IonLock™ Zinc increases singulation & adheres microbes to seed coats while providing valuable early emergence nutrition**
- **Builds stronger, healthier plants fast & massive roots**
- **Increases early-season iron, manganese and zinc uptake**
- **Multi-strain inoculant drives early-season root nodulation and N-fixation**
- **Puts soybeans in the position to manage season-long stressors and optimize genetic potential**
- **Increases biological activity in the soil**

Industry-Leading Bradyrhizobia Inoculant
Aggressive multi-strain inoculant to accelerate early nodulation


Revline
Proven, industry-leading biostimulant PPFM strains that generate massive root structures

Ether Enzyme Technology
Lipase & Mannanase enzymes that boost biological activity & activates nutrient availability

BioFertility & N-Fixing Microbes

Base Pail & Seed Fluency:

- 80/20 Talc/Graphite
- Mn & Fe
- 0.38 lbs IonLock Zinc equiv. to 1 qt 9% Zn EDTA



GUARANTEED ANALYSIS

BASE INGREDIENTS

Iron (Fe).....	0.70%
Manganese (Mn)	0.90%
Zinc (Zn).....	1.58%
1.58% Water Soluble Zinc	
Derived from Iron Oxide, Manganese Oxide and Zinc Sulfate.	

BIO-CAPSULE INGREDIENTS:

ALSO CONTAINS BENEFICIAL SUBSTANCES

<i>Azospirillum brasilense</i>	5.0x10 ⁷ CFU/g
<i>Azospirillum lipoferum</i>	5.0x10 ⁷ CFU/g
<i>Azotobacter chroococcum</i>	5.0x10 ⁷ CFU/g
<i>Azotobacter vinelandii</i>	5.0x10 ⁷ CFU/g
<i>Bacillus amyloliquefaciens</i>	5.3x10 ⁸ CFU/g
<i>Bacillus licheniformis</i>	1.1x10 ⁹ CFU/g
<i>Bacillus megaterium</i>	5.3x10 ⁸ CFU/g
<i>Bacillus pumilus</i>	8.8x10 ⁸ CFU/g
<i>Bacillus subtilis</i>	5.3x10 ⁸ CFU/g
<i>Paenibacillus azotofixans</i>	5.0x10 ⁷ CFU/g
<i>Thiobacillus ferrooxidans</i>	5.0x10 ⁷ CFU/g
<i>Trichoderma harzianum</i>	3.5x10 ⁵ CFU/g
<i>Bradyrhizobium elkanii</i>	1.00x10 ⁹ CFU/g
<i>Bradyrhizobium japonicum</i>	2.00x10 ⁵ CFU/g
<i>Methylobacterium hispanicum</i>	1.00x10 ⁹ CFU/g
Lipase.....	5.0x10 ⁴ μUnits/g
Mannanase.....	2.5x10 ⁹ μUnits/g

DIRECTIONS FOR USE & RECOMMENDED USE RATE

DO NOT ACTIVATE BIO-CAPSULE UNTIL READY FOR USE. Read and follow all label directions for use on package. Once ready to plant, remove the safety clips from the Bio-Capsules. Push down on the buttons atop the pail to release the biologicals into the base. Seal the pail and shake aggressively to blend the contents. Before each application, aggressively shake contents to ensure the correct volume is applied.

The combined products (1 pail) treat 40 units of soybeans at 140k seeds per unit.

PRODUCT DETAILS

Sold in 2-pail cases to treat 80 units of soybeans. Each pail is capable of treating 40 units of soybeans.